

PRODUCT INFORMATION

Anti-Spike (RBD) Antibody_AK3400

Description:

InVivo offers a series of different monoclonal antibodies for the detection of the SARS-CoV-2 Spike protein. All antibodies come from hybridoma cell lines that were generated by immunization with the recombinant full-length Spike protein (S) or the receptor-binding domain (RBD).

This antibody has been validated by ELISA and is specifically directed against an epitope that is located on the RBD. Additionally, ELISAs with directly coated antigens proved that this antibody also binds to the Alpha variant of the Spike protein and different RBD variants (Alpha, Alpha + E484K, Beta and Gamma).

This antibody was shown to neutralize the SARS-CoV-2 wildtype (Wuhan), as well as the Alpha and the Beta variant in cell culture. No other variants were tested.

All antibodies are produced exclusively under serum-free conditions from hybridoma and purified through one-step purification with Protein-A affinity chromatography.


Product-ID:	AK3400
Host:	Mouse
Clonality:	Monoclonal
Isotype:	IgG
Subclass:	mIgG2a/bk
Formulation:	Liquid, PBS, pH 7.4, 0.2 µm sterile filtered
Concentration:	≥ 0.5 mg/ mL
Purity:	≥ 90% (via analytical CGE under reducing conditions)
Conjugate:	Unconjugated
EC₅₀ RBD¹:	0.054 nM
EC₅₀ S1¹:	1.159 nM
EC₅₀ S¹:	0.374 nM
Specificity:	S (tested for: Wuhan and Alpha) RBD (tested for: Wuhan, Alpha, Alpha+E484K, Beta and Gamma)
Neutralization:	yes
IC₅₀ SARS-CoV-2 Wuhan²:	0.26 +/- 0.04 µg/ mL
IC₅₀ SARS-CoV-2 Alpha²:	0.12 +/- 0.01 µg/ mL
IC₅₀ SARS-CoV-2 Beta²:	0.12 +/- 0.002 µg/ mL

The product is for research use or for further manufacturing only.

¹ EC₅₀ values for 20 nM of coated antigen (Wuhan origin)

² Mean of two individual in-cell neutralization assays performed with Vera6 cells by our cooperation partner Dr. Valeria Falcone from the University Hospital Freiburg (Institute of Virology)


Purity (analytical CGE, under reducing conditions):




Peak	Size [kDa]	% of Total	Peak Identification
1	4.5	0.0	Lower Marker
2	5.7	0.0	System Peak
3	6.0	0.0	System Peak
4	7.9	0.0	System Peak
5	21.0	0.4	unknown
6	23.2	0.2	unknown
7	25.3	Peak Value 1	AK3400 Light Chain
8	52.9	0.9	unknown
9	57.9	Peak Value 2	AK3400 Heavy Chain
10	240.0	0.0	Upper Marker

Summation of peak values 1 and 2 results in a purity of $\geq 90\%$

Protein Activity (ELISA):






For ELISA analyses immunoassay plates were coated with different concentrations of SARS-CoV-2 S1 (10, 21, 42 and 84 nM). Antibody AK3400 was applied at concentrations ranging from 26.67 nM to 0.026 nM.

Bound antibody AK3400 was detected using Goat Anti-Mouse-IgG-HRP (10 ng/mL) and TMB. Determined EC₅₀ of antibody AK3400 for 20 nM S1 is 1.159 nM.

	84 nM S1	42 nM S1	21 nM S1	10 nM S1
EC ₅₀ [nM]	0.4551	0.5911	1.014	2.272





For ELISA analyses immunoassay plates were coated with different concentrations of SARS-CoV-2 Spike protein (3.6, 7, 14 and 29 nM). Antibody AK3400 was applied at concentrations ranging from 26.67 nM to 0.026 nM.

Bound antibody AK3400 was detected using Goat Anti-Mouse-IgG-HRP (10 ng/mL) and TMB. Determined EC₅₀ of antibody AK3400 for 20 nM Spike protein is 0.374 nM.

	29 nM S	14 nM S	7 nM S	3.6 nM S
EC ₅₀ [nM]	0.5491	0.5241	0.9229	1.992

Binding activity to SARS-CoV-2 Spike protein and receptor-binding domain variants (ELISA):



Neutralization Assay:

