



INVIVO

BioTech Services GmbH  
a BRUKER company



# Monoclonal Antibody Production

## APPLICATIONS:

- research use only
- Preclinical use
- Human *in vitro* diagnostics
- Veterinary diagnostics
- Food and feed safety

## INFORMATION:

- InVivo BioTech Services GmbH is certified according to
  - ISO 9001:2015
  - ISO 13485:2016
- Bulk production up to 1 kg

**InVivo BioTech Services GmbH is a leading European biomanufacturer with ISO 9001 and ISO 13485 certification. We have been a Bruker company since 2017.**

**For more than 25 years, InVivo has been committed to supporting academics and industries with our competence, expertise, and reliable products.**

**We are there with you to make you feel confident and to bring more efficiency when converting your ideas into practice.**

*In vitro* production in super spinner systems or bioreactors (2 L to 100 L) through high-cell density fermentation, occurs exclusively under serum-free conditions. For this purpose we use our proprietary medium **ISF-1** that allows cost-effective production processes and high product yields.

€ 999



# Welcome Offer

## Small-scale pilot production of monoclonal antibodies from hybridoma lines

### INCLUDED SERVICES:

- Mycoplasma testing
- Adaption to serum-free medium and cultivation in suspension (1 L scale)
- One-step purification by affinity chromatography
- Quality control by CGE and analytical SEC

### TURNAROUND TIME:

- 8 - 10 weeks

### DELIVERABLES:

- 10 mg purified antibody (guaranteed)
- Certificate of Analysis


### EXTRAS:

- Lyophilization
- Antibody labelling
- Aliquoting



**InVivo has cultivated more than 3500 different hybridoma lines so far, producing antibodies from all kinds of species (rat, mouse, rabbit, human, etc.) and of various subclasses (IgGs, IgMs, IgAs).**

Antibodies ordered in larger quantities (> 1g) and on a regular basis, qualify for our special **Antibody-on-Demand service**. This reduces costs and saves storage capacities at your facility. Just define your demand for 24 months, choose a batch size (min. 100 mg per call-off order) and delivery dates (fixed dates or flexible interval, e.g. monthly, quarterly, etc.) and we take care of the rest.



IVS CS-0151 Rev. 01

### InVivo BioTech Services GmbH

Neuendorfstraße 24a  
16761 Hennigsdorf  
Germany

Phone +49 (0) 3302 866 93-21/-22  
Fax +49 (0) 3302 866 93-62

[info.invivo@bruker.com](mailto:info.invivo@bruker.com)

For detailed product  
information please see  
[www.invivo.de](http://www.invivo.de)

