

This offer is valid for lyophilization of monoclonal antibodies and recombinant proteins.

Included Services:	<ul style="list-style-type: none">- Filling of lyophilization vials- Lyophilization- Labeling- Optical evaluation of the pellet
Deliverables:	<ul style="list-style-type: none">- Lyophilized product- Certificate of Analysis for liquid product before lyophilization

INSTRUCTIONS

Please complete this form and send it to info.invivo@bruker.com. Fields marked with an asterisk are mandatory. Not available or confidential information can be marked with "n/a".

CONTACT INFORMATION

	Billing Address	Delivery Address <i>(if different)</i>
Name*		
Company or Institution*		
Department		
Address*		
Phone*		
Email*		
VAT Number*		

ANTIBODY OR PROTEIN INFORMATION

Clone Name or Protein Name	
InVivo ID	
Total amount to be lyophilized	
Protein amount per Vial	mg / vial

Generally, the protein or antibody is concentrated to 1.00 – 2.00 mg/ mL before lyophilization and present in PBS Buffer pH 7.4. The 2 mL plastic tubes are usually filled up with 0.5 mL of the protein solution, the 4 mL glass vials with up to 1 mL and the 11 mL glass vials with up to 2 mL. Please let us know if you prefer other specifications.

ADDITIVES AND STABILISERS

Additives / Stabilisers	BSA	% (w/v)
	Sucrose	% (w/v)
	Other	% (w/v)
	Mixture of Additives:	

LYOPHILIZATION PROTOCOL

Protocol	InVivo Standard Other (Please provide protocol in attachment)
-----------------	------------------------------------------------------------------

LABEL

Clonename AK-Batch_LYO ___mg in ___ Date of Lyophilization YYYY-MM-DD 2...8°C InVivo BioTech Services GmbH

If special requirements for the label apply, then please provide information below.

VIALS

Choice of Vials	Glass vials with thread (I - IV) Plastic tubes of PP material with screw cap (V)
------------------------	-----------------------------------------------------------------------------------------




Figure 1: Choice of vials for lyophilization

PRICE FOR LYOPHILIZATION SERVICE (protein material not included):

695.00 € for up to 240 vials – applicable for all vials

1180.00 € for up to 735 vials – applicable for 4 ml glass vials and plastic tubes

1925.00 € for up to 1500 vials – applicable only for plastic tubes

STORAGE AND DELIVERY

Storage	Standard at 2-8 °C Other
Delivery	Standard at 2-8 °C Other